

CHALLENGER Series

SPECIFICATIONS

LOW WIDE

16 lbs*
21" Seat Width
19" Seat Height
26" Width x 28" Length
33" - 37" Handle Height
400 lbs Weight Capacity

WIDE

16 lbs*
21" Seat Width
21" Seat Height
26" Width x 28" Length
33" - 37" Handle Height
400 lbs Weight Capacity



THE CHALLENGER SERIES IS ONE OF OUR MOST VERSATILE WALKERS.

Small and compact, yet strong walker designed especially for versatility. Both models have 8" wheels and a comfortable, flexible back support.

FEATURES

ALUMINUM FRAME Superb stability and durability, yet lightweight

SEAT Solid plastic base with removable fitted waterproof cushion

WHEELS Solid grip and long lasting polyurethane tires, suitable for indoor & outdoor use

REAR FORK BRAKING SYSTEM

FLEXIBLE BACK STRAP
For added comfort

COLORS Vibrant Red, Black



OPTIONAL ACCESSORIES

- Cane Holder
- Slow Down Brakes
- Deluxe Back Strap Pad
- Oxygen Tank Holder
- Tote Bag
- Seat Mounted Tray
- Curb Climber
- Anatomical Hand Grips

*Net Weight of walker without accessories, backstrap and folding basket