

XPRESSO LITE Series

THE XPRESSO LITE
ALLOWS FOR A
MORE ENJOYABLE,
INDEPENDENT
LIFESTYLE.

Designed for the user that requires a light-weight frame, the Xpresso Lite series features low maintenance cable brakes, a soft, center folding padded seat and featherweight tires and wheels capacity of up to 350 lbs.



FEATURES

INNOVATIVE DESIGN Roman arch provides support, stability & style

SEAT Center folding durable waterproof covered foam cushion

COLLAPSIBLE BASKET
Foldable & removable

WHEELS 8" lightweight EVA tires

COLORS Ruby Red, Ebony Black



OPTIONAL ACCESSORIES

- Cup Holder
- Cane Holder
- Curb Climber
- Slow Down Brakes
- Deluxe Back Strap Pad
- Phone Holder
- Eyeglass Holder
- Mono Brake (Factory Installed Option)

SPECIFICATIONS

MINI

15 lbs*
18" Seat Width
18" Seat Height
22.5" Width x 26" Length
30.5" - 33.5" Handle Height
300 lbs Weight Capacity

REGULAR

15 lbs*
18" Seat Width
21" Seat Height
22.5" Width x 26" Length
30.5" - 34.5" Handle Height
300 lbs Weight Capacity

TALL

15 lbs*
18" Seat Width
24" Seat Height
22.5" Width x 26" Length
33.5" - 40" Handle Height
350 lbs Weight Capacity

*Net Weight of walker without accessories, backstrap and folding basket



Evolution
TECHNOLOGIES

Toll Free 1 800 556 2558

info@evolution.ca

www.evolution.ca