

XPRESSO WIDE Series



THE XPRESSO WIDE SERIES OFFERS THE MOST STABLE AND DURABLE PLATFORM FOR USERS WITH HEAVY-DUTY REQUIREMENTS.

Made with aerospace-grade aluminum and high strength steel, the Xpresso Wide series is designed with a weight capacity of 450 lbs.

FEATURES

INNOVATIVE DESIGN Roman arch provides support & stability

SEAT Center-folding durable waterproof covered foam cushion

BRAKE SYSTEM Cable-free/ Easy to adjust

COLLAPSIBLE BASKET Foldable & removable

WHEELS 8in heavy-duty polyurethane tires

FRONT & REAR FORKS CNC-milled aluminum forks for added strength

NEW PATENTED BACK STRAP Added comfort without compromising visibility

COLORS Red, Silver, Black



OPTIONAL ACCESSORIES

- Cup Holder
- Phone Holder
- Slow Down Brakes
- Eyeglass Holder
- Cane Holder

SPECIFICATIONS

WIDE

31 lbs*
20" Seat Width
22" Seat Height
24.8" Width x 25.5" Length
32.5" - 37.5" Handle Height
450 lbs Weight Capacity

LOW WIDE

30 lbs*
20" Seat Width
19" Seat Height
24.8" Width x 25.5" Length
30.5" - 35.5" Handle Height
450 lbs Weight Capacity

*Net Weight of walker without accessories, backstrap and folding basket



Toll Free 1 800 556 2558

info@evolution.ca

www.evolution.ca