

# XPRESSO ZERO Series



THE XPRESSO ZERO SERIES USES EVOLUTION'S PATENTED ROMAN-ARCH DESIGN FOR SUPERIOR STRENGTH AND STABILITY.

With the new padded seat cushion and wide back strap pad, users can sit in comfort and walk in style. The Xpresso Zero comes with an all-new braking system, allowing for effortless handle height adjustment and easy maintenance.

## FEATURES

**SEAT** 1.5" thick padded seat cushion for extra comfort

**BRAKE SYSTEM** Concealed cable / Easy to adjust and maintain

**COLLAPSIBLE BASKET** Fully covered for increased privacy. Foldable & removable

**WHEELS** 8" Lightweight EVA tires / Offset rear forks for wider gait

**WIDE BACK STRAP** Extra wide for added comfort

**HAND GRIPS** Removable design for easy maintenance

**COLORS** Carbon Black, Shiraz Red, Utopia Blue, Aquamarine Blue, Chrome Silver



## SPECIFICATIONS

### MINI

15 lbs\*  
17.5" Seat Width  
18" Seat Height  
25" Width x 23.5" Length  
28.5" – 31.5" Handle Height  
300 lbs Weight Capacity

### REGULAR

15 lbs\*  
17.5" Seat Width  
21" Seat Height  
26" Width x 23.5" Length  
31" – 35" Handle Height  
300 lbs Weight Capacity

### TALL

16 lbs\*  
17.5" Seat Width  
24" Seat Height  
27" Width x 24" Length  
33" – 37" Handle Height  
350 lbs Weight Capacity

## OPTIONAL ACCESSORIES

- Cup Holder
- Cane Holder
- Phone Holder
- Slow Down Brakes
- Eyeglass Holder
- Curb Climber
- Plastic Tray
- Oxygen Tank Holder
- Wire Basket

\*Net Weight of walker without accessories, backstrap and folding basket



Toll Free 1 800 556 2558

info@evolution.ca

www.evolution.ca