

EVOLUTION Series



SPECIFICATIONS

LOW

20 lbs*
17" Seat Width
19" Seat Height
21" Width x 30" Length
30" - 36" Handle Height
350 lbs Weight Capacity



REGULAR

23 lbs*
17" Seat Width
23" Seat Height
21" Width x 30" Length
30" - 38" Handle Height
350 lbs Weight Capacity



LOW WIDE

22 lbs*
21" Seat Width
21" Seat Height
26" Width x 30" Length
30" - 36" Handle Height
400 lbs Weight Capacity



WIDE

23 lbs*
21" Seat Width
23" Seat Height
26" Width x 30" Length
30" - 38" Handle Height
400 lbs Weight Capacity

AVAILABLE IN FOUR SIZES, FEATURES THE STRONGEST AND ONE OF THE MOST STABLE WALKERS ON THE MARKET TODAY.

The Evolution model uses 1-1/8" steel tubing, the same material used in mountain bike fabrication. Strong 1" tubing is used on the handles, every fastener is stainless steel or plated against corrosion. The Evolution series offers unmatched stability for users suffering from neurological conditions such as Multiple Sclerosis and Parkinson's or users requiring a high weight capacity.

FEATURES

STEEL FRAME Scratch resistant powder coated heavy-duty steel frame

SEAT Solid plastic base with removable fitted waterproof cushion

WIRE BASKET Large & removable

WHEELS 8" heavy-duty polyurethane tires

COLORS Black



OPTIONAL ACCESSORIES

- Cane Holder
- Curb Climber
- Slow Down Brakes
- Oxygen Tank Holder
- Mono Brakes
- In-Basket Tote Bag
- Seat Mounted Tray
- Forearm Platform System
- Arm Trough System
- Anatomical Hand Grips

*Net Weight of walker without accessories, backstrap and folding basket



Evolution
TECHNOLOGIES

Toll Free 1 800 556 2558

info@evolution.ca

www.evolution.ca