

# XPRESSO Series



THE XPRESSO, ALONG WITH ALL EVOLUTION WALKERS, WAS DESIGNED TO BE THE STRONGEST AND MOST STABLE OF ALL THE CENTER FOLDING WALKERS ON THE MARKET TODAY.

With cable-free brakes, a center folding seat, durable folding basket, the Xpresso folds easily and compact.

## FEATURES

**INNOVATIVE DESIGN** Roman arch provides support & stability

**SEAT** Center folding durable waterproof polyurethane cushion

**BRAKE SYSTEM** Cable-free / Easy to adjust

**COLLAPSIBLE BASKET** Foldable & removable

**WHEELS** 8" long lasting polyurethane tires, suitable for indoor & outdoor use

**COLORS** Ruby Red, Azure Blue, Ebony Black



## OPTIONAL ACCESSORIES

- Cup Holder
- Cane Holder
- Curb Climber
- Slow Down Brakes
- Deluxe Back Strap Pad
- Oxygen Tank Holder
- Wire Basket
- Eyeglass Holder
- Phone Holder

## SPECIFICATIONS

### MINI

18 lbs\*  
18" Seat Width  
18" Seat Height  
22.5" Width x 26" Length  
30.5" - 33.5" Handle Height  
350 lbs Weight Capacity

### REGULAR

18 lbs\*  
18" Seat Width  
21" Seat Height  
22.5" Width x 26" Length  
30.5" - 34.5" Handle Height  
350 lbs Weight Capacity

### TALL

18 lbs\*  
18" Seat Width  
24" Seat Height  
22.5" Width x 26" Length  
33.5" - 40" Handle Height  
350 lbs Weight Capacity

\*Net Weight of walker without accessories, backstrap and folding basket